

# Background Reduction Sheets

## Description:

Background Reduction Sheets obtains the best images of chemiluminescent Western blots and fluorescent blots and gels.

- **IMPROVED IMAGES** – Capture better images by removing background fluorescence and light noise from the environment
- **FLEXIBILITY** – Compatible with chemiluminescent and fluorescent Western blots, as well as gels stained with DNA and protein stains including SYBR green, SYBR Gold, ethidium bromide, SYPRO® Orange, SYPRO Ruby and more
- **INCREASED SENSITIVITY** – Improve signal to noise ratios
- **CONVENIENT** – Simply place the sheet under the blot or gel while imaging

## - Improve chemiluminescent Western blots

Plastic Wrap




Background Quenching Sheet




## - Improve chemiluminescent Western blots

Plastik Wrap



Background Reduction Sheets



## For Orders:

Catalog Number	Product	Size
NB057001010	Background Reduction Sheets	10 sheets