

# Performance Blotting Folders

## Description:

Performance Blotting Folders are transparent plastic folders that hold your blot while it is exposed to film or a CCD camera. These folders keep your blot flat for perfect film exposures, avoiding the wrinkling or bunching that can occur with plastic wrap.

The folders are not autofluorescent and are compatible with chemiluminescence and fluorescence imaging. The folders also protect blots during handling.

- **CLEAR IMAGES** – Avoid the wrinkling or bunching that can happen with plastic wrap
- **SUITABLE FOR ALL BLOTS** – Development folders are suitable for chemiluminescence and fluorescence imaging
- **SAFE BLOT HANDLING** – Move your blot without scratching or putting pressure on the surface
- **CONVENIENT SIZE** – Image up to 2 blots simultaneously in one folder

## For Orders:

Catalog Number	Product	Size
NB050702025	Performance Blotting Folders	25 folders
NB0507020A1	Performance Blotting Folders	100 folders