

Western Running Buffer, 10X

Description:

10X Western Running Buffer (10X) is containing Tris-glycine-SDS and the most commonly used buffer for sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE) of proteins.

- **Ready to dilute** – just use distilled deionized water for dilution
- **Reduced preparation time** – no reagents to weigh and filter
- **Quality-controlled reagent** – reproducible results

For Orders:

Catalog Number	Product	Size
NB001709100	Western Running Buffer, 10X	100 mL
NB001709500	Western Running Buffer, 10X	500 mL
NB001709102	Western Running Buffer, 10X	1 L